

Custom Component Tracking Solution

Introduction

ProcessHQ currently maintains an extensive online presence via its website – www.processHQ.com. The site contains information about the company, industry articles and links, and an online procurement solution

Custom Component Tracking

We have extended processHQ.com to specifically meet the growing need for managing custom components, such as hose assemblies. Each hose assembly is assigned a unique identifier (PHQ QC #) and recorded in our database. All component information, such as **manufacture date, customer purchase order number, part number**, etc. is correlated with and stored alongside the QC number. In addition, digital versions of pertinent documents, such as the quality assurance **material certificate, passivation certificate of conformance**, drawings, images, etc. are also viewable online. The database may be searched by the QC number or the customer PO number to locate the component of interest.

Thus, a precise record of all hose assemblies, including certificates, is always online and available for both viewing and printing. The component tracking database may be customized and integrated with your procurement solution to allow for single-click reordering of hose assemblies and other custom components. ProcessHQ offers complementary products in AdvantaLabel™ hose identification and HoseTrack™ RFID hose tracking and identification to further leverage the flexibility and power of our solution, facilitating an intelligent component replacement strategy.

High-Purity Component Procurement

The ProcessHQ web site allows you to easily navigate an extensive product database to quickly compare different products and select from a full line of specialized components, materials and equipment. The ProcessHQ online marketplace is specifically designed for the sanitary process manufacturing market with features such as custom hose assembly design and purchasing.

Security

The ProcessHQ site utilizes industry standard SSL encryption to insure privacy and security during network transmission. User names and passwords are never transmitted or stored in clear text. Furthermore, each hose assembly is assigned to a customer group, and only members of the group are granted access privileges to that assembly. Only the ProcessHQ site administrator is able to grant and revoke group membership, thus your company's components may only be accessed by authorized personnel as determined by your own policies and procedures.

Welcome Test Account from Berlex Laboratories to your custom components database.

Search For:

- pHQ QC #
- Your PO #

Search

List All

Refresh

Log out

QC # 100026	Part # 12-APSH-11-10-72	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 3/4 inch polyester braid reinforced silicone hose with Mini TC OE x 1.0 Inch TC other end, 72 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100037	Part # 16-APSH-10-10-72	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 1.0 inch polyester braid reinforced silicone hose with 1.0 inch TC each end, 72 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100039	Part # 16-APSH-10-10-120	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 1.0 inch polyester braid reinforced silicone hose with 1.0 inch TC each end, 120 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100048	Part # 16-APSH-10-10-144	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 1.0 inch polyester braid reinforced silicone hose with 1.0 inch TC each end, 144 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100073	Part # 12-APSH-11-11-72	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 3/4 inch polyester braid reinforced silicone hose with Mini TC each end, 72 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100084	Part # 12-APSH-11-10-72	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 3/4 inch polyester braid reinforced silicone hose with Mini TC OE x 1.0 Inch TC other end, 72 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100105	Part # 12-APSH-11-11-78	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 3/4 inch polyester braid reinforced silicone hose with Mini TC each end, 78 inch over all length		Customer Note:	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	
QC # 100167	Part # 32-APSW-N-10-10-16.5	Date Born: 1/23/03	PO # 8432764 OG	Add to PO
Description: 2.0 inch smooth wire reinforced silicone hose, 2.0 inch TC each end, 16.5 inch over all length(critical)		Customer Note: T-2901 DRN	More Info: AdvantaLabel (click to view text file) Material Certificate (click do download PDF file) Passivation Certificate (click do download PDF file)	