

Sanitary Hose Storage Rack

Protect your investment in sanitary flex hose assemblies with this versatile hose storage and organization solution. Designed specifically for bio-pharmaceutical clean room use, this robust, high-capacity unit is ideally suited for wet and corrosive environments. The rack facilitates drainability, organization, and safe transportation of a wide range of hose sizes and lengths.


Hoses of up to 5 feet in length are conveniently stored, sorted by size, in the hanging slots specifically sized for standard Tri-clamp ferrules. Longer hoses are placed on the top hanging rack for gentle and secure storage. The optional bottom shelf provides additional storage or functions as a drip tray. Customized designs are welcomed. All units are fabricated in the United States from high-grade 300-series stainless steel with sealed welds and a smooth bright finish, and feature sanitary, clean room casters for swivel and brakes for easy handling.

Key Features and Benefits

- Compact footprint, high capacity storage of over 60 flexible sanitary hoses
- Safe and secure handling and transportation of high-cost flex hose assets
- Ergonomic, easy access saves on labor
- Serves as an organization tool, sorting hoses by size or application
- Provides complete drainability
- Clean room quality fit and finish
- Sanitary casters (swivel / brakes)
- All 300-series stainless steel construction for corrosion resistance
- Optional utility tray and other customizations available
- Dimensions: 30" W x 30" L x 72" H

Storage Capacity

Hose Sizes	Capacity
1/2" - 3/4"	24
1" - 1-1/2"	24
2" - 2 1/2"	8
Top Rack	20



ProcessHQ
PURE AND SIMPLE